


NPGROLL



NEW POWER GENERATION

www.npgfitness.co.uk



**infrared
lights**



**profiled
trikes**



**16 items
of the massage**



**additional
accessorie**



intuitive panel B2




intuitive panel B2



	LXWXH	PRODUCT WEIGHT	MAX. USER WEIGHT	INFRARED POWER
J1	75 x 110 x 75 cm	93 kg	130 kg	300 W
B2	78 x 110 x 80 cm	75 kg	180 kg	500 W



- ✓ **breaking down fat tissue**
- ✓ **firming of the connective tissue**
- ✓ **improvement in blood supply**
- ✓ **release of muscle tension, fascia**
- ✓ **relief from excessive stress**
- ✓ **regression and prevention of orange peel skin**
- ✓ **stimulation and strengthening of muscles**



DEVICE PERSONALISATION

Upholstery colour (wide range selection available)



Logo embroidery



Trikes colours



White paint

Wooden

Black paint

Infrared Roll Massager

The operation is very simple – a specially profiled drum rotates, while rollers massage individual parts of the body, stimulating circulation and lymph flow. Infrared rays heat up the skin further improving circulation, metabolism and removing toxins from the body. Infrared rays strongly stimulate the immune system and accelerate collagen production.

Roll massage stimulates the lymphatic system to work, which results in the increased need of organs for oxygen and nutrients. stimulate the immune system and accelerate collagen production.

Roll massager is also a device for massaging the fascia. Excessive tension of the fascia leads to various disorders in the functioning of the whole organism. The fascia massage is perfect as a warm-up and complements training.



Roll massagers are ideal for:

- home,**
- solariums,**
- fitness clubs, gyms,**
- dedicated roll massaging parlours,**
- physiotherapy treatment rooms,**
- wellness and beauty salons,**



Sales:



Fitness Club 24
ul. Komorowicka 110
43-300 Bielsko-Biała
Poland

+48 33 486 90 07
sklep@fitnessclub24.pl

www.fitnessclub24.co.uk